

BIVENT MODE

Combination mode

Set 2 Pressure Levels (P High & PEEP)

Set 2 Pressure Supports (PS above P High & PS above PEEP)

Set 2 Time cycles (T High & T PEEP)

Patients can breathe spontaneously at any time during the time cycles

Indications for Bivent

Lung Recruitment

Pulmonary Contusions

ARDS

Patients who need inverse I:E ratio that you do not want to paralyze

****ALL changes on Bivent or any patient on inverse I:E ratio need to be discussed with the attending.